

Guanidinium thiocyanate (GITC)

lysis buffers for nucleic acid extraction



Lysis buffers to extract viral RNA

Guanidine thiocyanate is a potent chaotropic agent; thus, by interfering with the hydrogen bond network in aqueous solutions, it has a destabilising effect on macromolecules, especially proteins.

It is commonly used in cells and the lysis processes in virus particles to extract nucleic acids, as it denatures RNAse and DNAse enzymes¹ that would otherwise damage the extract.²

GITC lysis buffers to extract viral RNA are in growing demand, linked to the use of polymerase chain reaction (PCR) based assay.

Buffer composition (as reported by Scallan et al.3):

4 M Guanidinium thiocyanate (GITC)

55 mM* Tris-HCl

25 mM EDTA (Ethylenediaminetetraacetic acid)

3 % (v/v) Triton X-100

0.01 % (w/v) Bromophenol blue

(*NOTE: calculated from the total amount of 0.1 M Tris pH 7.6 added, diluted by the degree of volume expansion observed when the GITC goes into solution)

Method to produce one litre of 4M Guanidinium thiocyanate (GITC)/Triton X-100 Lysis buffer:

- 472.75 g of GITC is brought into solution initially by adding 400 ml of 0.1 M Tris HCl pH 7.6.
 This will require heating in a 65°C water bath and some shaking of the vessel (but with lid well secured). According to scientific literature, once fully dissolved the volume of the solution was 600 ml
- 2. Make up to 750 ml with 0.1 M Tris HCl pH 7.6
- 3. Add 50 ml of 0.5 M EDTA, mix
- 4. Add 30 ml Triton-X-100, mix

5. Volume made up to 1 L with 0.04 % (w/v) Bromophenol blue (DEPC-treated water can be used instead)

Note: Always read the chemical safety data sheet associated with the chemicals and carry out a full risk assessment.

Within our Thermo Fisher Scientific Laboratory Chemicals portfolio we offer all the necessary




Fisher Scientific Cat. No.	Description	Pack size
11352888	Bromophenol blue, ACS	1 g, 5 g, 25 g
11468716	Bromophenol blue sodium salt	10 g, 50 g
11308967	Bromophenol blue sodium salt, 0.04% w/v aq. soln.	100 mL, 500 mL
11329886	EDTA, disodium salt dihydrate,99+%, for molecular biology, DNAse, RNAse and protease free	100, 500 g, 2.5 kg
15450457	Ethylenediaminetetraacetic acid, electrophoresis grade, 99.4+%	100 g, 500 g, 2.5 kg
10503345	Guanidine thiocyanate powder assay: >=99.0 %	250 g, 1 kg
10741244	Guanidine thiocyanate, 99% (argentometric titration: >=98.5%)	100 g, 250 g, 1 kg
11454777	Guanidine thiocyanate, 99% assay (argentometric titration: \geq 98.5 to \leq 101.5%)	100 g, 250 g, 500 g
11975601	Guanidine thiocyanate, for molecular biology (argentometric Titration >=99.0 %)	100 g
15486719	Guanidine thiocyanate, molecular biology grade - assay (titration: 99.0% min.)	50 g, 250 g, 500 g
10790872	Hydrochloric acid, ACS reagent, ca. 37% solution in water	500 mL, 2.5 L
11463443	Hydrochloric Acid (Technical), Fisher Chemical	2.5 L
10163243	Tris Base, molecular biology grade	500 g, 1 kg, 5 kg, 10 kg, 25 kg
15446989	Tris(hydroxymethyl)aminomethane, ultrapure, 99.9%	1 kg, 5 kg
15855348	Tris, 99.0-101.0% (dry basis)	500 g, 1 kg, 5 kg, 10 kg, 25 kg
15845368	Tris, 99.8-100.1% (dry basis), molecular biology grade, ultrapure	500 g, 1 kg, 5 kg 10 kg
10671652	Triton™ X-100, 98%, for molecular biology, DNAse, RNAse and protease free	100 mL, 250 mL
15428059	Triton™ X-100, electrophoresis reagent	100 mL, 500 mL, 2.5 L
15456979	Water, endotoxin-free	500 mL, 5 L
10505854	Water, molecular biology grade	100 mL, 1 L, 4 L, 10 L, 20 L
10191231	Water, for molecular biology, DNAse, RNAse and Protease free	1 L, 5 L
15805408	Water, RNAse-free, DEPC treated, molecular biology grade, ultrapure	10 x 1 mL, 25 mL, 100 mL, 500 mL, 1 L

Bulk and semi-bulk quantities available

Thermo Fisher Scientific is an approved supplier for several pharmaceutical organizations globally. Several products, including **Guanidinium Thiocyanate** (product codes: **41111**, **B21250**, **43394**, **J65104**) are available in bulk and semi-bulk quantities in full compliance with local regulations.

For more information eu.fishersci.com/go/thermo-scientific-chemicals





Trust us to help manage your fine chemical purchasing process for:

- Bulk
- Semi Bulk
- Specialty Chemicals

Discover more and request a quote for specialty and semi-bulk quantities at **thermofisher.com/chemicals**

References

- 1. McGookin R (1985) RNA extraction by the guanidine thiocyanate procedure. Methods Mol Biol 2:113-6.
- 2. Chomczynski P, Sacchi N (1987) Single-step method of RNA isolation by acid guanidinium thiocyanate-phenol-chloroform extraction. *Anal Biochem* 162(1):156-9.
- Scallan M F, Dempsey C et al. (2020) Validation of a Lysis Buffer Containing 4 M Guanidinium Thiocyanate (GITC)/ Triton X-100 for Extraction of SARS-CoV-2 RNA for COVID-19 Testing: Comparison of Formulated Lysis Buffers Containing 4 to 6 M GITC, Roche External Lysis Buffer and Qiagen RTL Lysis Buffer. bioRxiv ePub: 1-6 (This article is a preprint and has not been certified by peer review)

For more information visit our website:

eu.fishersci.com/go/thermo-scientific-chemicals

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk



© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.