Instruction Manual Manuel d'instructions

AGITATEUR MAGNETIQUE CHAUFFANT

F20500413FI - F20500162FI

General Information / Informations Générales

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Caution hot surface!
Attention surface chaude!

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.

This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Fisher Scientific reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Safety Regulations / Regles de Securite

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Beware of the effect of the magnetic field on cardiac pacemakers and data media.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it. It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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Contents / Index

EN

1.	INTRODUCTION	3
2.	ASSEMBLY AND INSTALLATION.....	3
3.	START-UP	3
	3.1 USING THE THERMOREGULATOR VTF (ONLY FOR F20500413FI)	3
	3.2 ERROR MESSAGES	3
4.	END-OF-WORK OPERATIONS	3
5.	MAINTENANCE	4
	5.1 CLEANING.....	4
6.	TECHNICAL DATA	4
7.	ACCESSORIES	4

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The application of new technology has led to the creation of these modern magnetic stirrers for laboratory use whose basic features are safety, reproducibility of results, high performance and energy saving.

The structure of the unit is made out of pressure die-cast aluminium treated with epoxy resins offering a high resistance to the many chemical aggressions typical of the laboratory environment.

The meticulously studied design of the pressure die-cast casing combines aesthetic fluidity with ergonomics. A wide central recess on both sides offers a good grip and facilitates bench-top handling whilst the handle moulded directly into the rear of the die-cast structure can be used to transport the instrument once it has been disconnected from the power supply.

The control panel is easily accessible, is positioned at a safe distance from sources of heat and is protected against possible damage caused by liquid spills thanks to a special run-off groove.

The pressure die-cast structure is designed so that accidental liquid spills cannot reach the internal parts of the unit.

The electrical connections are recessed into the rear of the instrument and offer optimum electrical protection of the external electrics in the case of liquid spills.

The engineering of the unit facilitates internal inspection if technical assistance is required.

The heating plate is made out of aluminum alloy with special coating and is specially treated to ensure:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Optimum temperature homogeneity across the plate
- High resistance to thermal stress and thermal shock, chemical aggressions, scratches and surface abrasions.

NOTE: using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

The stirring speed is adjusted by an impulse microprocessor which makes it possible to obtain optimum speed progression. Speeds of from 50 to 1.500 rpm can be selected using the speed-selector knob.

The heating plate temperature is controlled by an analog bulb thermoregulator fitted directly onto the plate itself.

The heating plate temperature and all the set temperatures are thus controlled constantly.

2. Assembly and installation

The unit is delivered in special shock-resistant packaging and comes complete with: Power cable

Check that the values on the rating plate correspond to those of the power supply.

Connect the unit to the power supply using the power cable supplied.

Turn the unit on using the mains switch, the switch will light-up indicating that the unit is powered.

3. Start-up

Turn the unit on using the mains switch and place a beaker with magnetic stir bar on the heating plate. Select the stirring speed using the speed-selector knob. Stirring speeds of from 50 to 1500 rpm can be selected.

Select the heating plate temperature using the temperature control knob. Temperatures of between 50 and 370 °C can be selected.

3.1 Using the thermoregulator VTF (only for F20500413FI)

Screw the threaded support rod (optional accessory, A00001069) into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the receptacle making sure that it is completely immersed in the sample. Connect the two instruments (F20500413FI) by plugging the VTF into the dedicated socket on the back of the stirrer.

Select the operating temperature required on the VTF thermoregulator. Turn the temperature control knob on the front panel of the stirrer to maximum.

The heating magnetic stirrer always has primary control of the heating plate temperature.

When using the VTF thermoregulator always select the max temperature on the stirrer. The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the temperature setting on the stirrer meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting.

3.2 Error messages

The unit is fitted with safety devices which cut-off the power supply to the heating plate in the case of malfunctions.

The led display on the front panel indicates the type of malfunction:

Display

Flashing Stirrer LED (once every second)

Flashing Heating LED (once every second)

Flashing Heating LED (2 times every second)

Cause

Fault in the stirring system

Overheated heating plate or thermocouple open

Fault in the temperature reading circuit

Should any of the above occur, please contact your service centre.

4. End-of-work operations

When the unit is not in use it is advisable to turn off the mains switch, the light will go out.

5. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. It may be necessary to replace the fuses. The unit is fitted with two fuses located in the socket on the back of the unit. To replace one or more of the fuses, disconnect the unit from the power supply and use a screwdriver to lift the fuse box cover. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletised).

5.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent when the plate is cold.

6. Technical data

Power supply	230 V / 50-60 Hz o 115V / 60 Hz (+/-10%)
Dimensions WxHxD	165x115x280 mm
Weight	2.5 kg
Power input	630 W
Heat output	600 W
Diameter of the heating plate	155mm
Programmable temperature range	50 – 370 °C
Type of temperature control	Analog
Overtemperature protection	Yes
Motor rating output	10 W
Stirring capacity	15 liters of H ₂ O (F20500162FI) 20 liters of H ₂ O (F20500413FI)
Programmable speed range	100 – 1500 rpm
Counter-reaction	Constant speed even when the viscosity changes
Type of motor control	Electronic
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Pollution degree CEI EN61010-1	2
Overvoltage category	II
Noise level	<< 80 dBa
Max altitude	4000 m

7. Accessories

F208B0063FI	Electronic Vertex VTF -10...+300 °C
40000781*	Spiral electric cable complete of plug
A00001069	Support rod
A00001056	Magnetic stir bar diam. 6x35mm
A00001061	Magnetic stir bar diam. 9,5x60 mm
A00001060	Magnetic stir bar diam. 10x40mm
A00001071	Hemispheric bowl for 250 ml flasks
A00001072	Hemispheric bowl for 500 ml flasks
A00001073	Hemispheric bowl for 1000 ml flasks

* For connection to other contact thermometers other than VTF

		
F208B0063FI	40000781	A00001069
		
A00001056	A00001061	A00001060
		
A00001071	A00001072	A00001073

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