## ***Instruction Manual*** ***Manuel d'instructions***

### **VTF Fuzzy Logic Thermoregulator**

#### **General Information / Informations Générales**

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#### **This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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#### **This unit has been designed and manufactured in compliance with the following standards:**

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**UL 61010-1**

**CAN/CSA-C22.2 No.61010-1**

Fisher Scientific reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The Thermoregulator is ideal for many applications and meets the most demanding requirements in terms of precision, reliability and flexibility of use thanks to the application of Fuzzy Logic technology. The Fuzzy Logic electronic automatically adapts thermoregulation to the varying factors such as power, load and thermal dispersion specific to each application, by optimizing both overheating and the oscillations around the temperature set-point.

The accuracy and precision of thermoregulation at each end the scale and whatever the volume being processed, is a fundamental characteristic of the Vertex.

**WARNING:** when the Thermoregulator is connected to an instrument other than the VELP Scientifica AREX, AREC.X or PW10, please check the compatibility of the DIN 5 pole connector as illustrated in wiring diagram in this manual.

**NOTE:** the most precise results are obtained when measuring the temperature of aqueous solutions.

The microprocessor offers various other functions:

- a working time of up to 24 hours and 59 minutes can be set with automatic switch off;
- the maximum sample temperatures reached during the test can be recorded.

The structure is made of a non-scratch technopolymer resistant to chemical agents and offers a high level of IP54 electrical protection in compliance with regulation CEI EN 60529.

The instrument has an in-built safety and control circuit which is constantly active and shuts-down the thermoregulator immediately in the following situations: the temperature probe is not connected; the temperature probe is faulty (cut-off or short-circuited) and/or the temperature is out-of-range.


Fig. 1

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Thermoregulator complete with: power cable suitable for connection to VELP Scientifica magnetic stirrers AREX and AREC.X and to the VELP Scientifica PW10

- Pt100 temperature probe
- Instruction manual

### 2.1 Connecting the Pt100 probe

Connect the probe to the Thermoregulator as shown on Fig. 2.

### 2.2 Installing the VTF

The VTF has an innovative integrated system to simplify installation on the support rod and facilitate the positioning of the Pt100 probe in the most commonly used containers. The two clamps allow height-regulation whilst a sliding clamp allows horizontal regulation (Fig. 3).


Fig. 2


Fig. 3

### 2.3 Using the VTF with magnetic stirrers

The power cable is suitable for connection to the most common heating magnetic stirrers with a dedicated socket for the remote control of the heating plate such as the VELP Scientifica models AREX (code F20520163) and AREC.X (code F20500061). The socket must supply a tension of between 9 and 16V DC. The integrated installation system allows the instrument to be installed on support rods with a diameter of between 10 and 13 mm.

### 2.4 Other applications

When using the thermoregulator for other applications bear in mind that the transistor output can take a maximum current of 50mA. Where power loads of up to 2200W are required the use of an external power relay is necessary (code A00000001, VELP model PW10). For the thermoregulation of liquids that are not compatible with the construction material of the probe (stainless steel AISI 316), the use of a glass-coated probe is recommended (code A00000003).

In cases in which it is necessary to thermoregulate a liquid that is not in close proximity to the VTF, a 1m probe extension cable is available (code A00000002).

### 3. Operating controls

Install the thermoregulator on the support rod and place the probe in the liquid to be processed. To power the VTF plug the spiral cable into the dedicated socket on the heating magnetic stirrer and turn the stirrer on. The display shows the software version after which it flashes for approximately five seconds showing the last temperature Set Point used (Fig. 4).

Fig. 4

The current Pt100 probe temperature reading is then displayed on the main window.


Fig. 5

Thermoregulation starts automatically based on the last temperature Set Point used (Fig. 5).

**IMPORTANT:** the magnetic stirrer always exercises primary control of the heating plate temperature. To control the heating plate temperature using the VTF thermoregulator, set the magnetic stirrer temperature to maximum. The magnetic stirrer temperature control can also be used as a safety thermostat since the heating plate will not exceed the temperature set on the stirrer. In this case a longer heating time will be necessary in order for the sample to reach the set point temperature of the VTF.

#### 3.1 Setting the working temperature

From the main window, press the central  button. The set-up window is displayed (Fig. 4).

Set Point Temperature                      From -10 to +300 °C                      Default value: 40 °C

With the display flashing use the  and  keys to select the temperature required. The temperature setting is saved when no keys are touched for approx. 5 seconds (Fig. 5).

If no key is pressed for approx. 5 seconds when the set-up window is flashing, the value shown on the display is automatically saved. After 5 seconds the VTF will evaluate the temperature reading of the probe immersed in the liquid and will proceed to thermoregulate the liquid to the temperature selected.

**NOTE:** to display the maximum temperature reached since the instrument was turned on, with the main window displayed keep the  key pressed. The maximum temperature reached is deleted from the memory every time the instrument is turned off.

#### 3.2 Setting the timer

From the main window press the  key twice to display the working time.

Working time                      From 00:00 to 24:59 (h:min.)                      Default value: 0:00 (Fig. 6)


Fig. 6

With the display flashing use the  and  keys to select the working time required.

To switch from hours to minutes press the  key once or simply wait 4 seconds (Fig. 7).

Fig. 7

With the display flashing use the  and  keys to select the working time required. The value is saved when no keys are touched for approx. 5 seconds. If no key is pressed for 5 seconds the display returns the main window and the previously displayed value is saved.

Count-down starts the moment the working time has been saved. The alternate flashing of the two arrows to the left of the display indicates that the VTF is running in timer-mode (Fig. 8).


Fig. 8

At the end of the working time the instrument automatically turns off the thermoregulation, and the display shows "End" (Fig. 9).


Fig. 9

To return to the main window in order to start new work-cycle, press . The working time will return to the default value of 0:00 and the maximum temperature reached can be viewed by pressing . The working time can be modified when the VTF is in use even if the instrument is running in timer-mode. If the working time is set to 0:00 the instrument will run in continuous mode.

**NOTE:** to display the time left, with the main window displayed keep the  key pressed.

The time left is automatically reset to zero every time the instrument is turned off.


### 3.3 Temperature alignment

The micro-chip carries out automatic probe calibration therefore a probe calibration procedure is unnecessary. Nevertheless, an offset value of from - 9.9 to + 9.9 °C can be set.

Temperature alignment                      From - 9.9 to + 9.9 °C                      Default value    0.0 °C

Turn the instrument on by keeping the  and  keys pressed.

Use the  and  keys to modify the “Temperature alignment” parameters (Fig. 10).


If no key is pressed for 5 seconds the display shows the software version and the last set point value.

If no key is pressed for a further 5 seconds the display shows the probe temperature reading and saves the value.

### 3.4 Unit of measure

Turn the instrument on by keeping the  and  keys pressed.

Use the  and  keys to modify the “Unit of measure” (Fig.11).

If no key is pressed for 5 seconds the display shows the software version and the last set point value. If no key is pressed for a further 5 seconds the display shows the probe temperature reading in the new unit of measure and saves the value.

### 3.5 Safety devices and alarm signals

The instrument has an in-built safety and control circuit which is constantly active and shuts-down the VTF when an alarm message appears on the display:

- AL .1 = Unconnected probe or faulty probe (cut-off or short circuited) or temperature out of range.
- AL .2 = Slow temperature increase read by the probe. This function is active in case of:
  - Probe temperature < 50°C
  - T set point – T probe > 5°C
- AL .3 = Fast temperature decrease read by the probe.

To reset the thermoregulator turn it off and on again after having found and removed the cause of the alarm.

In order to bypass the AL .2 and AL .3, turn the instrument on by keeping the  and  keys pressed and select AL.oF using  and  keys. Press  to confirm.

## 4. Start-up

- Install the VTF on the support rod so that the probe is immersed in the liquid by at least 15mm making sure that it does not make contact with the bottom of the container;
- Make sure the magnetic stirrer is turned on;
- Power the VTF by plugging the spiral cable into the dedicated socket on the heating magnetic stirrer;
- Set the magnetic stirrer temperature to maximum;
- Select the temperature and the working time required;
- To optimize thermoregulation and improve temperature homogeneity, stir the sample gently but continuously during operation.

### WARNING

- When working with set points near the temperature of evaporation, check that the probe remains immersed in the sample by at least 10mm throughout the work-cycle.
- The sample being processed is subject to evaporation and for this reason it may not reach the temperature selected.

## 5. End-of-work operations

To interrupt thermoregulation and leave the VTF thermoregulator powered, turn off heating plate using the knob on the magnetic stirrer. In this way heating stops but the VTF continues to display the temperature.

If the timer has been set, the thermoregulator will automatically stop at the end of the set time and the display will show “End”.

To reactivate thermoregulation press .

At the end of the process, if the VTF has to be left connected to the magnetic stirrer, turn the stirrer off.

## 6. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

### 6.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent.

## 7. Technical data

### GENERAL

		Default
Power supply		12 V (9÷15)
Power	W	1.2
Absorption	mA	100 max
Output	mA	50 max
Size	mm (LxHxD)	75x145x120
Weight	Kg	0.250
Display		4 digit LCD
Type of thermoregulation		Multi function Fuzzy Logic
Support rod	∅ mm	From 10 to 13
Connector		5 pole 270° DIN
Construction material		Polypropylene (PP)
Temperature range	°C	0...+50
Storage temperature range	°C	- 10...+ 60
Max. humidity	%	90%

### FUNCTIONS

Temperature settings	°C / °F	- 10...+ 300 / 14...572
Selection interval	°C / °F	1 / 1
Resolution	°C / °F	0.2 / 1
Precision	°C / °F	± 0.5 / ± 0.9
Working time settings	h:min	From 0:00 to 24:59
Continuous mode		Possible
Time left reading		Possible
Maximum temperature reached reading		Possible

### PROBE


Type of probe		Pt100
Probe speed τ90%	sec	5
Probe length	mm	250
Immersion depth	mm	15 minimum

### SAFETY

Probe short-circuit	Error message and shut-down	automatic
Break in probe circuit	Error message and shut-down	automatic
Out of temperature range	Error message and shut-down	automatic
Level of electrical protection	IP54	CEI EN60529

## 8. Accessories


A00000001	Derivation element PW10
A00000002	Probe extension cable, length 1m
A00000003	Glass probe
A00000280	Clamp with probe support

## 9. Spare Parts

40000579	Pt100 probe	40000580	Friction assembled
----------	-------------	----------	--------------------

## 10. Wiring Diagram

## 11. Declaration of conformity

We, the manufacturer Fisher Scientific, under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1/2001  
EN 61326-1/2006  
2011/65/UE (RoHS)  
2012/19/UE (RAEE)

and satisfies the essential requirements of the following directives:

2006/42/EC  
2006/95/EC  
2004/108/EC

plus modifications and that the documents listed in the annex I are available at Velp's offices as per the machinery directive

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