



---

Important Notes..... 6

- 1. General Practice ..... 6
- 2. Safety ..... 6
- 3. Maintenance of the instrument ..... 7
- 4. Declaration of Conformity ..... 7

Chapter I Introduction ..... 8

Chapter II Features..... 9

- 1. Normal working conditions ..... 9
- 2. Transportation and storage conditions..... 9
- 3. Basic parameters ..... 10
- Table 1: 100-120V..... 10
- Table 2: 200-2400V..... 11

Chapter III Preparation Work ..... 12

- 1. Structure diagram ..... 12
- 2. Display panel diagram ..... 12
- 3. Instructions for buttons ..... 13

Chapter IV Operation Guidance..... 14

- 1. Examination before startup ..... 14
- 2. Startup..... 14
- 3. Setting temperature and heating time ..... 14
- 4. Start and stop instrument..... 15
- 5. Temperature measurement and calibration..... 15

Chapter V Failure Analysis and Troubleshooting ..... 16

Chapter VI Accessories..... 17

Warranty..... 18

---

**This manual covers the models shown below**

| <b>NA Model</b> | <b>EU Model</b> | <b>Voltage</b> | <b>Description</b>               |
|-----------------|-----------------|----------------|----------------------------------|
| 88860021        | N/A             | 100-120V       | FS Drybath Stdrd 1 blck 100-120V |
| 88860022        | N/A             | 100-120V       | FS Drybath Stdrd 2 blck 100-120V |
| 88860023        | N/A             | 100-120V       | FS Drybath Stdrd 4 blck 100-120V |
| 88860024        | 15387928        | 200-240V       | FS Drybath Stdrd 1 blck 200-240V |
| 88860025        | 15397928        | 200-240V       | FS Drybath Stdrd 2 blck 200-240V |
| 88860026        | 15317938        | 200-240V       | FS Drybath Stdrd 4 blck 200-240V |



**IMPORTANT!** Before using this product, read this entire operation manual carefully. Users should follow all of the operational guidelines contained in this manual and take all necessary safety precautions while using this product. Failure to follow these guidelines could result in potentially irreparable bodily harm and/or property damage.

Caution all internal adjustments and maintenance must be performed by qualified service personnel.

Material in this manual is for information purposes only. Fisher Scientific is committed to a continuing program of product development and improvement, and reserves the right to change information, such as specifications, appearance, and dimensions, described in this document without notice. Fisher Scientific makes no representations or warranties with respect to this manual. In no event shall Fisher Scientific be held liable for any damages, direct or incidental, arising out of or related to the use of this manual.

No part of this manual may be reproduced or transmitted in any form or by any means, including photocopying, recording, or using information storage and retrieval systems, for any purpose other than the purchaser's own use, without the express written permission of the manufacturer.

Any other product names and services identified in this manual are trademarks or registered trademarks of their respective owners. No such use, or the use of any trade name, is intended to convey endorsement or other affiliation with Fisher Scientific.

©2017 Fisher Scientific. All rights reserved.

---

All rights reserved. Without prior written consent from Thermo Fisher Scientific, any part of the manual should not be duplicated, reproduced or translated into any other languages.

**Caution:**

All rights reserved. Thermo Fisher Scientific reserves the right to modify this manual at any time without notice.

**Before the initial use of this instrument, please carefully read this manual.**

---

# Important Notes

## 1. General Practice

### Caution:

Important information is contained in any item and should be carefully read. Failure to observe the instruction would result in damage or abnormal function of the instrument.

### Warning!

The warning message requires extremely careful operation of a certain step. Failure to observe the instruction may result in serious personal injury.

## 2. Safety

During operation, maintenance and repair of this instrument, the following basic safety notes should be observed. In case of failure to follow these instructions, the warnings or notes indicated herein, the basic protection provided by the instrument, its safety criteria of design and manufacture, and its predicted use range would be impaired. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

### Caution:

This instrument is for indoor use.

### a) The ground connection

To avoid the electric shock, the input power line must be reliably grounded. The instrument is equipped with the three-pin plug that has the third pin (the pin connects the ground), therefore, the plug should be used with the grounded power socket only. This is a kind of safety device. If the plug cannot be plugged into the power socket, please ask the electrician to install a correct power socket, so as to make the grounded plug work for safety.

Thermo Fisher Scientific is not responsible for any injury as a result of the user's failure to observe the following requirements.

### b) Servicing and Replacement of Parts

The operator should not take the instrument apart without permission. Only qualified professionals are allowed to replace elements or adjust components inside the instrument. Replacement of components while unit is connected to power and/or turned on is prohibited.

### c) Observe voltage/power requirements

Before the instrument is connected to the AC power source, the voltage of the power source should be the same with the required voltage of the instrument (a deviation of  $\pm 10\%$  is allowed). The rated load of the power socket should not be lower than the requirements of the instrument.

### d) The power cord

The instrument should use the power line cord attached to it. If the power line is destroyed, it must be replaced but not be repaired. The replacement should be carried out with the power line of the same type and same specifications. No items should be put on the power cord when the instrument is in operation.

### e) Connecting and disconnecting to/from power source

The user should hold the plug to remove from power source. When connecting the plug, user should make sure it has been fully plugged in; when removing the plug, do not pull the power line forcefully.

### f) Placement of the instrument

This instrument should be fixed in a low RH and low dust place away from water source (e.g. sink or water pipeline) and the room should be well ventilated, and free of corrosive gas or interference of strong magnetic field. The instrument should not be placed in wet or location. The openings on this instrument are for ventilation circulation and in order to avoid over-heat of the instrument body, they should not be blocked or covered. When a single set of instrument is used, the interval

between ventilation opening before and after the instrument and its nearest object should not be less than 25cm. Also, don't use the instrument on loose or soft surface, or the air inlet of its bottom might be blocked. Excessive temperature will impair performance or result in failure of the instrument. This instrument should not be used in locations subjected to direct sun light.

The instrument should be kept away from hot gas, oven and all other heat sources. If the instrument is to be stored for a long time, the power plug should be withdrawn and the instrument covered with soft cloth or plastic film to avoid entrance of dust. The product is powered by connecting the mains plug to a standard socket-outlet. Always place the product in such a way that it is easy for the operator to disconnect the product from the mains supply.

#### g) Explanation of symbols



Attention, read user manual before use.

#### Warning!

There is a sign of "CAUTION: HOT SURFACE! or ATTENTION: SURFACE HOT" on the instrument. The metal part (module) near the sign should not be touched with any part of the body when the Instrument is operating in a high temperature state or just finished operation to avoid burns!



Conductor Terminal



Alternating Current Protective

#### Caution:

In any of the following cases, immediately cut off the power supply, withdraw the power plug from the power socket and contact the supplier for service:

- Liquid drops into the inside of the instrument.
- The instrument is rained on or water is spilled on it.
- The instrument works abnormally, especially when generating an abnormal sound or odor.
- The instrument is dropped or its casing is damaged.
- The function of instrument obviously changes.

### 3. Maintenance of the instrument

The holes of the blocks should be regularly cleaned with the damp cloth to ensure the test tube be well contacted with the wall of the holes so as to have good heat conduction. If the surface of the instrument is polluted, it can be cleaned with a slightly damp soft cloth.

#### Warning!

When cleaning the instrument, the power supply should be shut off and unplugged. The instrument surface should be cleaned with a non corrosive cleaning agent.

### 4. Declaration of Conformity

We hereby declare under our sole responsibility that this product conforms to the technical requirements of the following standards:

|            |                             |
|------------|-----------------------------|
| CEEMC:     | EN 61326-1                  |
| CE Safety: | EN 61010-1 CE               |
| Safety:    | EN 61010-2-010              |
| UL:        | 61010-1/CSAC22.2 NO.61010-1 |
| ROHS2.0:   | 2011/65/EU                  |

---

# Chapter I Introduction

The product is a dry bath with advanced microcomputer control, which can be widely applied to sample reservation, enzyme reservation and reaction, DNA amplification, electrophoresis degeneration and serum coagulation, etc. The characteristics of the product are as follows:

- The digital display for control of the temperature.
- Heat blocks help avoid possible contamination while containing the sample inside a tube.
- The heat blocks are easy to replace, clean and disinfect and are suitable for various tubes sizes.
- Built-in over-temperature protection device warrants sample and user safety.
- Temperature deviation adjustment.

---

# Chapter II Features

This chapter introduces the usage, transportation, storage conditions of the instrument, as well as its basic parameters, performance and functions.

## 1. Normal working conditions

Ambient temperature: 10°C ~ 30°C

Relative humidity: ≤80%

Power supply: 100-120V~ 50/60 Hz or  
200-240V~ 50/60 Hz

Operation: Up to 2000m

Pollution degree: 2

Overvoltage category: II Indoor use

### Caution:

Before using the instrument, please make sure the working condition meets the above requirements. Especially that the power line is reliably grounded.

## 2. Transportation and storage conditions

Ambient temperature: -20°C ~ +55°C

Relative humidity: ≤80%

### 3. Basic parameters

Table 1: 100-120V

|                         |   |                          |                                  |
|-------------------------|---|--------------------------|----------------------------------|
| NA Catalog Number       | 88860021  | 88860022                 | 88860023                         |
| Model Name              | Digital 1 block 100-120V                                    | Digital 2 block 100-120V | Digital 4 block 100-120V         |
| Temperature Range       | Ambient Temperature +5°C ~ 130°C (Ambient Temperature 25°C) |                          |                                  |
| Temperature Uniformity  | ≤ ±1.0°C  |                          |                                  |
| Temperature Accuracy    | ≤ ±0.5°C@37°C   |                          |                                  |
| Temperature Fluctuation | ≤ ±0.5°C  |                          |                                  |
| Temperature Rise Time   | ≤20 min (rise from 30°C to 130°C)                           |                          | ≤25min (rise from 30°C to 130°C) |
| Size (LxWxH)            | 288x200x100 mm  | 318x200x100 mm           | 450x200x100 mm                   |
|                         | 11.3x7.87x3.94 in   | 12.5x7.87x3.94 in        | 17.7x7.87x3.94 in                |
| Weight                  | 2.9 kg 6.39 lbs.  | 3.3 kg 7.28 lbs.         | 4.7 kg 10.36 lbs.                |
| Fuse Protector (Φ5x20)  | F 250V 2.5A   | F 250V 5A                | F 250V 8A                        |

**Table 2: 200-2400V**

|                         |  |                          |                                    |
|-------------------------|--|--------------------------|------------------------------------|
| NA Catalog Number       | 88860024   | 88860025                 | 88860026                           |
| EU Catalog Number       | 15387928   | 15397928                 | 15317938                           |
| Model Name              | Digital 1 block 200-240V                                   | Digital 2 block 200-240V | Digital 4 block 200-240V           |
| Temperature Range       | Ambient Temperature +5°C ~ 13°C (Ambient Temperature 25°C) |                          |                                    |
| Temperature Uniformity  | ≤ ±1.0°C   |                          |                                    |
| Temperature Accuracy    | ≤ ±0.5°C@37°C  |                          |                                    |
| Temperature Fluctuation | ≤ ±0.5°C   |                          |                                    |
| Temperature Rise Time   | ≤20 min (rise from 30°C to 130°C)                          |                          | ≤ 25 min (rise from 30°C to 130°C) |
| Size (LxWxH)            | 288x200x100 mm   | 318x200x100 mm           | 450x200x100 mm                     |
|                         | 11.3x7.87x3.94 in  | 12.5x7.87x3.94 in        | 17.7x7.87x3.94in                   |
| Weight                  | 2.9 kg 6.39 lbs.   | 3.3 kg 7.28 lbs.         | 4.7 kg 10.36 lbs.                  |
| Fuse Protector (Φ5x20)  | F 250V 2.5A  | F 250V 2.5A              | F 250V 3.15A                       |

**Caution:** The instrument only use for heating operation: from low temperature to high temperature.

# Chapter III Preparation Work


This chapter introduces structure of the instrument, user interface and functions of all buttons and preparation before startup. Read the content in the chapter carefully before using the instrument for the first time.

## 1. Structure diagram



- 1. Housing
- 2. Front Cover
- 3. Display Window
- 4. Product Label
- 5. Power Socket, Fuse Housing, Power Switch
- 6. The Handle, to lift hot blocks

## 2. Display panel diagram



---

### 3. Instructions for buttons



Setting or adjusting button: Press the button to display present setting temperature and press the button again to display present time. After the unit has reached set point and temperature is stable, calibration can be set by pressing the button.



Under the status of setting, press "▲" or "▼" to modify the flashing figure.



The instrument will enter temperature calibration mode if press the two buttons simultaneously under the status of stopping.



When setting or adjusting, pressing the button can move the flashing figure. During the operation, the display can be switched between the temperature and the time.



For starting up or stopping program.

---

# Chapter IV Operation Guidance

This chapter introduces the displaying window and the operation of buttons of the dry bath with constant temperature.

## Notes:

The illumination pattern of the 'Heat Light' indicates the Operation status.

- Heat Light OFF (Unit Stopped) = "Set/Cal" button used for Setpoint adjustment.
- Heat Light blinks/ flashes as Block is heating or cooling toward setpoint.
- Heat Light ON (Unit Started and stable at setpoint) = "Set/Cal" used for Calibration.

## 1. Examination before startup

Please confirm the following content before plugging the plug into the power socket.

1. The power source is in accordance with the voltage requirement of the instrument (please refer to the chapter II for the requirement of the voltage);
2. Make sure the plug has been fully plugged into the power socket;
3. The power line is reliably grounded.

## Caution:

If the display of the instrument is abnormal after startup, please turn off the power source immediately and contact with the supplier.

## 2. Startup

Press the power switch and the instrument will power on. All indicator lights and display serial ports will be turned on. About 3 ~ 5 seconds later, temperature display window will show the real-time temperature. All indicator lights will be turned off.

## 3. Setting temperature and heating time

Set Temperature and Hold-Time Setpoint:

1. Connect power cord to appropriate electrical supply.
2. Make Power switch, located near power cord on rear, to the 'on' position.
3. Confirm the Heat-indicator is 'off'. If 'on' or 'blinking', press "Stop".
4. Set-up Temperature:
  - A. Press the "SET/CAL" button. The instrument is in the Temperature Setpoint mode with the display window showing the previous temperature set value, and the left digit flashing.
  - B. Press "▲" to increase or "▼" to decrease the flashing digit of the Setpoint.
  - C. Press "SHIFT" to move the flashing right, to the next digit.
  - D. Repeat Steps 4B and 4C to complete Temperature setpoint selection.
5. Setting Hold-Time Setpoint:
  - A. After the temperature is selected, press "SET/CAL" button. The instrument is in the Hold-time Setpoint mode with the display window showing the previous Hold-time set value, and the left digit flashing.
  - B. Press "▲" to increase or "▼" to decrease the flashing digit of the Hold-Time. Choose a time from 00:01 to 99:59, hours and minutes.
  - C. Press "SHIFT" to move the flashing right, to the next digit.
  - D. Repeat Steps 5B and 5C to complete Hold-Time setpoint selection. For Continuous operation, Choose 99:59 then increase to display "- :- -".
  - E. After the Time setting is complete, Wait for 5 seconds and the instrument will automatically exit the Setpoint mode. The display will now show the current measured temperature.

6. Press "Start" button to begin heating the Block. The green "Heat Light" will blink as Block is heating toward set point.
7. Allow 30-minutes for temperature to stabilize. The heat indicator stops blinking and remains constantly ON to indicate the temperature is stable, and instrument is ready to use.

**Instruction:**

There is limitation for temperature and time setting. The setting temperature range is: 15 ~ 135°C; the time range is : 0~99h59min and continuous mode "- :- -".

## 4. Start and stop instrument

**START:** Press "START/STOP" to start operation, the indicator light "HEAT" will be flashing. When the instrument achieves the setting temperature, the "HEAT" light will stop flashing. This indicates that the unit will enter hold temperature status. The digital display window will display the real-time temperature and the timing function by pressing "SHIFT" button. When the holding time equals to the setting time, the instrument will alarm 5 beeps.

**STOP:** Under the working status, press "START/STOP" to stop working. At that time, corresponding working indicator light "HEAT" will be turned off. The digital display window displays the real-time temperature.

## 5. Temperature measurement and calibration

1. The instrument has been operational for minimum of 30-minutes at desired setpoint.
2. An independent temperature measurement device is in the designed 5mm hole in Block, along with paraffin oil, or similar heat-transfer oil.
3. Confirm the Heat-indicator is 'ON' steady, not 'blinking'.

4. Compare the Display and independent measurements. And check if the temperature gap is >0.2°C or not.
5. If less than 0.2°C, no calibration should be made.
6. If more than 0.2°C, calibrate the Temperature by these steps:
  - A. Press and hold the "SET/CAL" button. The instrument is in the Temperature Calibration mode with the display window showing the measured temperature value, and the left digit flashing.
  - B. Press "▲" to increase, or "▼" to decrease, the flashing digit to match the same digit of the independent measurement.
  - C. Press "SHIFT" to move the flashing right, to the next digit.
  - D. Repeat Steps B and C to complete Temperature calibration.
  - E. After the temperature calibration is completed, press "SET/CAL" button to exit the calibration mode and continue the Start mode.
7. The Heat-indicator may start 'blinking' as the temperature is stabilized to the calibrated temperature.
8. When the Heat-indicator is 'ON' steady, the temperature is stable.
9. After the instrument stable operated for additional 10 minutes, compare the display temperature and measurement again. If the temperature gap is >0.2°C, repeat all steps of 6.

**Caution:**

1. To ensure the accuracy of temperature, the instrument shall be calibrated after temperature stabilizes for 30 minutes.
2. Calibrate the instrument with qualified Grade II standard filled thermometer.
3. Calibration point: The middle hole of the heat block. Fill with paraffin oil and soak the thermometer bulb.

---

# Chapter V Failure Analysis and Troubleshooting

This chapter introduces possible failures of this instrument, explanation and troubleshooting.

**Caution:**

The user is not allowed to open the housing of the instrument for inspection during the warranty period. If any failures that require opening the housing for inspection occur, the supplier or manufacturer should be contacted timely.

## Failure Analysis and Troubleshooting

| NO. | Error  | Cause                                | Soultion  |
|-----|--|--------------------------------------|---|
| 1   | Display window is not turned on after switching on power supply                  | Power supply fails to be connected   | Check power supply and connect it                                 |
|     |  | Fuse burned out                      | Replace fuse  |
|     |  | Damage of the switch                 | Replace the switch  |
|     |  | Others                               | Contact Thermo Fisher Scientific                                  |
| 2   | Temperature display window shows "001" and the buzzer alarm in a di...di...sound | Short circuit of sensor              | Check if the connecting line of sensor is damaged, and replace it |
| 3   | Temperature display window shows "002" and the buzzer alarm in a di...di...sound | Open circuit of sensor               | Check if the connecting line of sensor is damaged, and replace it |
| 4   | Temperature display does not conform to actual temperature                       | Sensor is damaged or in poor contact | Contact Thermo Fisher Scientific                                  |

# Chapter VI Accessories

## Accessories

| NA Cat. No. | EU Cat. No. | Heating Blocks   | No. of Holes |
|-------------|-------------|--|--------------|
| 88860101    | 15327938    | For 6mm test tubes   | 46           |
| 88860102    | 15337938    | For 10mm test tubes  | 28           |
| 88860103    | 15347938    | For 1.5ml test tubes   | 28           |
| 88860104    | 15357938    | For 2.0ml test tubes   | 28           |
| 88860105    | 15367938    | For 12mm & 13mm test tubes   | 24           |
| 88860106    | 15377938    | For 15mm & 16mm test tubes   | 15           |
| 88860107    | 15387938    | For 17mm & 18mm test tubes   | 12           |
| 88860108    | 15397938    | For 20mm test tubes  | 8            |
| 88860109    | 15307948    | For 25mm test tubes  | 6            |
| 88860110    | 15317948    | For 0.5ml test tubes   | 40           |
| 88860111    | 15327948    | For 96 well Elisa plate  | 1            |
| 88860112    | 15337948    | For mixed size test tubes 32 holes for 6mm test tubes<br>21 holes for 10mm test tubes                            | 32+21        |
| 88860113    | 15347948    | For mixed size test tubes 18 holes for 1.5ml test tubes<br>10 holes for 2.0ml test tubes                         | 18+10        |
| 88860114    | 15357948    | For mixed size test tubes 3 holes for 25mm test tubes<br>12 holes for 13mm test tubes 6 holes for 6mm test tubes | 3+12+6       |
| 88860115    | 15377948    | For mixed size test tubes 30 holes for 0.5ml test tubes<br>20 holes for 0.2ml test tubes                         | 30+20        |
| 88860116    | 15387948    | For 15 ml flat bottom test tubes   | 15           |
| 88860117    | 15397948    | For 50 ml flat bottom test tubes   | 4            |
| 88860118    | 15307958    | For 15ml conical bottom test tubes   | 15           |
| 88860119    | 15317958    | For 50ml conical bottom test tubes   | 4            |
| 88860120    | 15337958    | For 96 well none skirted PCR plate   | 1            |
| 88860121    | 15347958    | For 96 well half/full skirted PCR plate  | 1            |

\* Dry bath blocks can be autoclaved.

---

# Warranty

When used in laboratory conditions and according to these operation instructions and maintenance, this product is warranted for 24 months against defective materials or workmanship. The 24 month warranty period begins from the delivery date of this product.

For product quality or performance issues, contact Fisher Scientific Customer Service.

## North America

United States  
1-800-766-7000  
fishersci.com

Canada  
1-800-234-7437  
fishersci.ca

## Europe

Austria:  
+43(0)800-20 88 40  
at.fishersci.com

Netherlands:  
+31(0)20 4887 70 00  
nl.fishersci.com

Belgium:  
+32(0)56 260 260  
be.fishersci.com

Norway:  
+47 22 95 59 59  
fishersci.no

Denmark:  
+45 70 27 99 20  
fishersci.dk

Portugal:  
+351 21 425 33 50  
pt.fishersci.com

Germany:  
+49(0)180 5258221  
de.fishersci.com

Spain:  
+34 002 239 303  
es.fishersci.com

Ireland:  
+959(0)1 885 5854  
ie.fishersci.com

Sweden:  
+46 31- 68 94 30  
fishersci.se

Italy:  
+39 02 950 59 478  
it.fishersci.com

Switzerland:  
+41(0)56 618 41 11  
ch.fishersci.com

Finland:  
+358(0)9 8027 6280  
fishersci.fi

UK:  
+44(1)1509 555 500  
fisher.co.uk

France:  
+33(0)388 67 14 14  
fishersci.fr

Find out more at **fishersci.com**

© 2017 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.