Fisherbrand
QUALITY. RELIABILITY. VALUE.

Instruction Manual - Gebrauchsanweisung - Manuel d'utilisation - Manual de funcionamiento

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Table of content

Bottle Top Dispenser

• Intended use of the instrument	5
• Specifications	5
• General safety instruction	5
• Diagram of the Bottle Top Dispenser	5
• Restriction of use	5
• Operating instructions	5
• Autoclaving	5
• Troubleshooting Guide	5
• List of recommended reagents	5
• Chemical Resistance Chart at 20°C	5
Intended use of the instrument	6
Specifications	6
General safety instruction	6
• A) Dispenser when not in use :	6
• B) Dispenser when in use :	6
Diagram of the Bottle Top Dispenser	7
Restriction of use	8
• Before using the Dispenser	8
• Assembly	8
Operating instructions	8
• Priming	8
• Dispensing	8
• Volume adjustment knob	9
• User calibration procedure	9
• Maintenance / Cleaning	9
• Dis-assembling the Dispenser for Cleaning and Servicing	10
Autoclaving	13
• Dis-assembling for Autoclaving	13
• Re-assembling after autoclaving	13
Troubleshooting Guide	14
List of recommended reagents	15
Chemical Resistance Chart at 20°C	16

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General safety instruction

A) Dispenser when not in use :

- Piston Barrel should always be empty
- Nozzle cap should be fixed on the nozzle

B) Dispenser when in use :

- Nozzle cap should be removed
- Place a receiving vessel under nozzle before starting the operation
- Never use force

Diagram of the Bottle Top Dispenser

Restriction of use

NEVER use the Dispenser with :

- Liquids which are not compatible with PTFE, FEP and Borosilicate Glass
- Hydrofluoric acid
- Liquids which contain solid particles
- Temperature limits are 15°C to 40°C

Before using the Dispenser

Check that the instrument has not been damaged in transit

Assembly

The Dispenser is packed with the dispense nozzle attached and the inlet feed tube removed. The length of FEP inlet tubing provided should be adjusted to fit the particular reservoir. Longer lengths on inlet tube are available on request. The threaded platform base of the Dispenser has a 30mm screw thread. The assembled Dispenser is screwed to the reservoir using gentle hand torque applied to the threaded platform base only. Removal should also be by means of hand torque applied to this same base. Do not operate the piston until the unit is safely and fully mounted on the reservoir bottle. Six adaptators are supplied to suit containers with a 28mm, 30mm, 32mm, 36mm, 40mm and 45mm screw neck.

Operating instructions

Priming

Place a container under the Dispenser's delivery nozzle. Remove the Nozzle Cap. Set the "Volume adjustment knob" to the maximum volume to enable free movement of the piston.

Prime the unit with a few gentle up and down strokes, taking the piston right down to its lowest stop position and lifting it up. Repeat until a steady bubble free flow is visible in the barrel.


Dispensing

- Ensure that the nozzle cap is removed.
- Ensure that a receiving container is placed the nozzle.
- Dispenser is now ready for dispensing.

Volume adjustment knob

Volume adjustment knob is simple and easy to use. It is screw type knob. Simply unscrew the knob and slide up or down to set the desired volume by aligning the pointer on the volume adjustment knob with the scale on the main housing.

Tighten the screw after aligning with the scale.

Volume adjustment knob

User calibration procedure


Dispenser has been laboratory calibrated to its nominal volume. However, due to changes in environmental conditions and the viscosity of the media which you dispense, a re-calibration might be required. You can either re-calibrate at regular intervals such as once a week or whenever you notice that the dispensed volume is different from the volume displayed by the unit.

To fully re-calibrate your Dispenser follow the following steps:

1. Set the Dispenser to the nominal volume or any other volume which is the most common volume you dispense.
2. Follow the common rules for calibration used in statistical quality control (ISO 8655/2). Set the volume and Dispense five full volumes of distilled water at 20°C on Electronic Balance to establish the actual mean volume of liquid dispensed. If the gravitational average result varies from the volume displayed, you should re-calibrate the Dispenser.
3. For re-calibration pull the cap outwards to expose the Calibration Nut.
4. Using the calibration tool, turn the calibration nut clockwise to reduce the volume and anticlockwise to increase the volume. Repeat this procedure a few times till the desired volume is achieved.

Cap

Maintenance / Cleaning

Note: All maintenance should be carried out wearing suitable eye protection and protective clothing. If in doubt, consult your safety officer.

1. Make sure that the Dispenser is completely empty.
2. Place the instrument into an empty sink together with its reservoir.
3. Unscrew the threaded platform base from the reservoir and lift the dispenser's intake tube carefully out of the reservoir, whilst tapping it against the reservoir's aperture to shake off any droplets from

the intake tube.

4. Hold the dispense nozzle over the aperture of the reservoir and apply gentle piston strokes in order to return any syringe contents into the reservoir.
5. Empty the instrument completely and flush thoroughly with distilled water.
6. If the piston barrel is still not completely clean, you need to dis-assemble the dispenser. Refer Dis-assembling procedure given below.



Dis-assembling the Dispenser for Cleaning and Servicing

Procedure to dis-assemble the PISTON

- 1- Pull the cap outwards to expose the Calibration Nut.
- 2- Unscrew the Calibration Nut with the help of calibration tool to dis-assemble the Piston and shaft out of the main housing.

- 3- Unscrew the piston from the shaft

Procedure to dis-assemble the BARREL

- 1- Unscrew the Barrel Cap using the special tool provided with the instrument.

- 2- Remove the Barrel Cap.

- 3- Lift the upper housing to expose the Barrel and cover.

- 4- Unscrew Chuck Nut 1 and remove Barrel cover.

- 5- Glass Barrel is now exposed.

- 6- Gently pull the barrel upwards to detach it from the Valve Manifold.
- 7- Barrel has been dis-assembled for cleaning.

Procedure to dis-assemble the DELIVERY and the VALVE MANIFOLD

- 1- Remove the delivery pipe cover by pulling it upwards from the slot.

- 2- Unscrew and remove the delivery pipe.

- 3- Unscrew the chuck Nut 2.


- 4- Remove Chuck Nut 2 and pull out the valve manifold

Autoclaving

Dis-assembling for Autoclaving

- 1- Unscrew Chuck Nut 1

- 2- Pull Chuck Nut 1 along with barrel cover, upper housing and piston all the way up.

- 3- Autoclave the two sub-assemblies at 0 121 C and 15 psi pressure for 10-15 mins.

Re-assembling after autoclaving

- 1- Push the Piston in the Glass Barrel gently and go all the way down. (Caution : Ensure that the Nozzle cap is removed)

- 2- Tighten chuck nut 1 by screwing it properly.

- 3- Dispenser is now ready for use. No Re-calibration is required after Autoclaving. However, a quick calibration check is recommended.

Troubleshooting Guide

Most problems can be avoided by reading and following the instructions in this operating manual. Below we list some of those most commonly experienced along with suggestions for solving them. If, however, these should not resolve the issue, or if you have questions not covered below, please contact Fisher Scientific.

Trouble	Possible cause	Solution
Air bubbles appear in discharge nozzle	<ul style="list-style-type: none"> Liquid reservoir is empty 	<ul style="list-style-type: none"> Refill reservoir and prime unit
	<ul style="list-style-type: none"> Too fast filling action 	<ul style="list-style-type: none"> Fill and dispense more slowly
	<ul style="list-style-type: none"> Leaking piston 	<ul style="list-style-type: none"> Clean piston. If problem persists replace piston
	<ul style="list-style-type: none"> Leaking discharge valve 	<ul style="list-style-type: none"> Clean unit by flushing thoroughly. If problem persists, replace the discharge valve
Barrel does not fill with liquid	<ul style="list-style-type: none"> Inlet tube not fitted correctly 	<ul style="list-style-type: none"> Connect inlet tube correctly
Dispensing not possible	<ul style="list-style-type: none"> Blocked dispense nozzle 	<ul style="list-style-type: none"> Disassemble the dispense nozzle and flush with cleaning fluid
	<ul style="list-style-type: none"> Discharge valve stuck 	<ul style="list-style-type: none"> Clean unit by immersing valve assembly in cleaning fluid. If problem persists, replace valve assembly
Wrong dispense volume	<ul style="list-style-type: none"> Instrument not calibrated 	<ul style="list-style-type: none"> Follow steps of User Calibration

List of recommended reagents

1,4-Dioxane	Dichloroethane	Sulphuric acid, 98%
1-Butanol	Lactic acid	Diethylene glycol
Acetaldehyde	Magnesium chloride	Dimethylformamide (DMF)
Acetic acid, 50%	Mercury Chloride	Ethanol
Acetone	Methanol	Ethy acetate
Acetonitrile	Methyl propyl ketone	Ethylene glycol
Acrylonitrile	Monochloroacetic acid	Formaldehyde, 40%
Adipic acid	n-Amyl acetate	Formic acid, 100%
Allyl alcohol	n-Butyl acetate	Glycerol
Aluminium chloride	Nitric acid, 70%	Heating oil (Diesel oil)
Amino acids	Nitrobenzene	Hexane
Ammonium chloride	Octane	Hydrochloric acid, 37%
Ammonium hydroxide, 30%	Oleic acid	Iodine / potassium iodide solution
Amyl alcohol (Pentanol)	Oxalic acid	Isobutanol
Aniline	Pentane (n-)	Isopropanol (2-Propanol)
Barium chloride	Perchloric acid, 10%	Isopropyl benzene (Cumene)
Benzaldehyde	Phenol	Tartaric acid
Benzene (Benzol)	Phosphoric acid, 85%	Tetrachloroethylene
Benzine (Gasoline)	Potassium chloride	Tetrahydrofuran (THF)
Benzyl alcohol	Potassium dichromate	Tetramethylamminium hydroxide
Biuret reagent	Potassium hydroxide	Toluene
Boric acid, 10%	Potassium permanganate	Trichloroacetic acid
Calcium carbonate	Propanol	Trichloromethane (Chloroform)
Calcium chloride	Propylene glycol	Triethylene glycol
Carbon tetrachloride	Propylene Oxide	Turpentine
Chromic acid, 10%	Salicylic acid	Urea
Chromic acid, 50%	Silver acetate	Xylene
Chromosulfuric acid	Silver nitrate	Zinc chloride, 10%
Copper sulfate	Sodium acetate	Zinc sulfate, 10%
Cresol	Sodium dichromate	
Dichlorobenzene	Sodium hydroxide, 30%	

These recommendations are carefully checked and correspond to the current state of knowledge. If you need statements for chemicals which are not given in the list, please do not hesitate to contact us.

CAUTION

- Do not use HF or reagents not compatible with PTFE or Borosilicate Glass.
- If used with strong acids, it is advised to rinse & remove instrument at the end of every working day & store safely.

Chemical Resistance Chart at 20°C

Liquids dispensed with the dispenser will be in contact, constantly, with the following materials; Borosilicate glass, (BSG), PTFE & FEP. The following table is a guide to help with the queries regarding liquid compatibility.

Please note that these tables are just a guide. We recommend that if there is a question regarding liquid compatibility you should exercise caution in use and refer to other chemical tables where available. Good laboratory practice would be to rinse out the liquid handling unit at the end of each day with distilled water to prevent corrosive liquids being left in contact with the parts for too long.

CHEMICAL Acids	BSG	PTFE	FEP
Acetic, Glacial	R		
Acetic, 25%	R	R	R
Hydrochloric, Concentrated	R		
Hydrochloric, 25%	R	R	R
Sulphuric, concentrated	R		
Sulphuric, 25%	R	R	R
Nitric, Concentrated	R		
Nitric, 25%	R		
Phosphoric, 25%	R	R	R
Formic, 25%	R	R	R
Trichloroacetic 10%	R	R	R
Formic Acid, 85%	R	R	R
Arsenic Acid	R		
Boric Acid, 10%	R	R	R
Chromic Acid, 20%	R	R	R
Hydrofluoric Acid, 35%	NR	exceptions	R
Phosphoric Acid 85%	R	R	R
Nitric Acid, 50%	R	R	R
Sulphuric Acid, 95%	R	R	R
Alkalies			
Ammonium Hydroxide, 25%	R	R	R
Sodium Hydroxide	R	R	R
Potassium Hydroxide	R	R	R
Sodium Hydroxide	R	R	R
Alcohols			
Methanol, 98%	R	R	
Ethanol, 98%	R		
Ethanol, 70%	R		

CHEMICAL Acids	BSG	PTFE	FEP
Isopropanol, n-Propanol	R		
Amyl Alcohol, Butanol	R		
Benzyl Alcohol	R	R	R
Ethylene Glycol	R	R	R
Propylene Glycol	R	R	R
Glycerol	R	R	R
Hydrocarbons			
Hexane, Xylene	R	R	R
Toluene, Benzene	R	R	R
Kerosene, Gasoline	R		
Tetralin, Decalin	R		
Halogenated Hydrocarbons			
Methyl Chloride	R		
Chloroform	R	R	R
Trichloroethylene	R	R	R
Monochlorobenzene, Freon	R		
Carbon Tetrachloride	R	R	R
Ketones			
Acetone	R	R	R
Methyl Ethyl Ketone	R	R	
Isopropylacetone	R		
Methyl Isobutyl Ketone	R		
Ethyl Acetate	R	R	
Methyl Acetate	R		
Amyl & Propyl Acetate	R		
Butyl Acetate	R	R	R
Propylene Glycol Acetate	R		
2-Ethoxyethyl Acetate	R		
Methyl Cellosolve Acetate	R		
l Benzoate	R		
Isopropyl Myristate	R		
Tricesyl Phosphate	R		

	CHEMICAL Acids	BSG	PTFE	FEP
Oxides-Ethers				
	Ethyl Ether	R		
	1,4 Dioxane & Tetrahydrofuran	R	R	R
	Dimethylsulphoxide (DMSO)	R	R	R
	Isopropyl Ether	R		
Solvents with Nitrogen				
	Dimethyl Formamide	R	R	R
	Diethylacetamide	R	R	
	Triethanolamine	R		
	Aniline	R	R	R
	Pyridine	R	R	R
Miscellaneous				
	Phenol, Aqueous, 10%	R		
	Formaldehyde Solution, 30%	R	R	R
	Hydrogen Peroxide, 30%	R	R	R
	Silicone Oil & Mineral Oil	R		
	Pyridine	R	R	R
	Acetaldehyde	R	R	R
	Ammonia, 25% aq. Sol	R	R	
	Ammonium	R		
	Calcium Chloride aq. Sol	R	R	R
	Chlorine	R	R	R
	Chlorobenzene	R		
	Fluorinated Hydrocarbons	R		
	Hexane	R	R	R
	Iodine (tincture of)	R	R	
	Potassium Chloride aq. Sol	R		
	Potassium Chloride aq. Sol	R		
	Potassium Permanganate aq. Sol	R		
	Magnesium Chloride aq. Sol	R		
	Methylene Chloride	R	R	R
	Sodium Carbonate	R		
	Sodium Dichromate	R	R	R

CHEMICAL Acids	BSG	PTFE	FEP
Phenol, 100%	R	R	R
Mercury	R	R	R
Silver Nitrate	R	R	R
Toluene	R	R	R
Hydrogen Peroxide, 30%	R	R	R
Xylene	R	R	R
Zinc Chloride, 10%	R	R	R
Zinc Sulphate, 10%	R	R	R

Key

R = resistant **VR** = virtually resistant **SR** = slightly resistant **NR** = non-resistant

Exceptions = resistant with exceptions

Notes : depends on temperature - up to 300°C

Fisherbrand[®]
QUALITY. RELIABILITY. VALUE.

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Flaschenaufsatz-Dispenser

12867913
2.5 mL

12877913
5.0 mL

12887913
10 mL

12897913
30 mL

12807923
50 mL

12817923
100 mL

Flaschenaufsatz-Dispenser

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The Bottle Top Dispenser is a general purpose laboratory instrument intended for use in laboratories for dispensing reagents and chemicals which are compatible (Please refer the compatability chart on page XX) with the instruments.

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Vol.-Bereich	Skala	Genauigkeit		CV	
				±%	±mL
±% ±ml	VK	0.6	0.015	0.2	0.005
±% ±ml	0.1	0.5	0.025	0.1	0.005
0,25–2,5 ml	0,05 ml	0,6	0,015	0,2	0,005
0,5–5 ml	0,1 ml	0,5	0,025	0,1	0,005
1–10 ml	0,2 ml	0,5	0,050	0,1	0,010
2,5–30 ml	0,5 ml	0,5	0,150	0,1	0,030
5–60 ml	1,0 ml	0,5	0,300	0,1	0,060
10–100 ml	2,0 ml	0,5	0,500	0,1	0,100

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