Instruction Manual

Fisher Scientific Compact Overhead Stirrer

English
page 3 to 13


Français
page 15 to 25

Español
page 27 to 37

Deutsch
page 39 to 49

ENGLISH	3
<i>Safety instructions</i>	4
<i>General instructions</i>	5
<i>Approvals</i>	6
<i>Packing list</i>	6
<i>Description</i>	7
<i>Technical data</i>	8
<i>Installation in the laboratory</i>	9
<i>Environmental operating conditions</i>	9
<i>General operating instructions</i>	10
<i>Routine cleaning, service, transport & storage</i>	12
<i>Warranty and product liability</i>	13

Safety Instructions

Each operator should read all instruction before using this device.

The stirrer must be securely fixed to a stable support. Mount to a Fisher Scientific stand (14-500-355 / 1198-4985) using a Fisher Scientific heavy duty clamp (14-500-357 / 1290-7106). If other stands or clamps are used, the unit must be secured such that it will not fall if tipped.

CAUTION: This is not an explosion proof stirrer. Do not use in a highly flammable or explosive environment.

Spinning paddles or impellers can cause severe personal injuries. Operators must use extreme care and good judgment when mixing at any speed.

All mixing paddles and impellers must be in good condition with straight shafts. If the stirrer vibrates at high speeds check the paddle shaft for damage and repair or replace if needed.

Extreme care must be taken when mixing chemicals to ensure that no chemicals are splashed outside the mixing vessel. Care must be taken when changing to faster mixing speeds.

Ensure that the mixing impeller does not contact the containment vessel.

Do not operate while standing in water. Keep the unit dry and do not immerse any part, except the mixing paddle into any liquids. Protect from splashing.

Ensure that no loose clothing, jewelry, or hair can become entangled in any rotating parts. A fast spinning chuck can cause injury to operator.

Power can be interrupted to the stirrer by pressing the power button. In an emergency, disconnect the mains cord.

General Instructions

Read safety instructions.

Unpack the stirrer carefully. Ensure the stirrer is not damaged and verify the contents against the packing list. If the stirrer is damaged or the contents do not correspond with the packing list, notify your supplier immediately.

Read the manual carefully! Make sure that all users read this manual thoroughly before operating the stirrer. For further questions or inquiries, please contact:

Fisher Scientific

United States

Customer Service: 1-800-766-7000

To Fax an Order: 1-800-926-1166

To Order online: www.fishersci.com

Canada

Customer Service: 1-800-234-7437

To Fax an Order: 1-800-463-2996

To Order online: www.fishersci.ca

Europe

Austria: +43(0)800-20 88 40

Belgium: +32 (0)56 260 260

Denmark: +45 70 27 99 20

Germany: +49 (0)180 5258221

Ireland: +353 (0)1 885 5854

Italy: +39 02 950 59 478

Finland: +358 (0)9 8027 6280

France: +33 (0)3 88 67 14 14

Netherlands: +31 (0)20 487 70 00

Norway: +47 22 95 59 59

Portugal: +351 21 425 33 50

Spain: +34 902 239 303

Sweden: +46 31-68 94 30

Switzerland: +41 (0)56 618 41 11

UK: +44 (0)1509 555 500

Approvals

North America

Conforms to ANSI/UL Std No. 61010-1 2nd Ed.
Certified to CAN/CSA-22.2 Std No. 61010-1 2nd Ed.
& 61010-2-051-04

Europe

EN61326 (EN55011, conducted and radiated
emissions, EN61000-3-2, EN61000-4-2, EN61000-4-3,
EN61000-4-4, EN61000-4-5, EN61000-4-6)

CSA	Canadian Standards Association
ANSI	American National Standards Institute
ISA	Instrument Society of America
CE	Communauté Européenne (European Community)
ITS	Intertek Testing Services (www.worldlab.com)
ETL	Electrical Testing Laboratory


Packing List

Compare the contents of the shipment to the items on the following list to ensure all parts are received with the stirrer. Do not discard the container and packing materials until all parts are accounted for. It is recommended that the original packing materials be kept for transportation and storage.


- Compact Overhead Stirrer with chuck Part # (NA 14-500-211 / 1520-3807)
- Power Supply with 115 volt cord set or (2) 220 volt cord sets
- Support Bar (including allen screw, lock washer and key)
- Chuck Key
- Instruction Manual

Compact Overhead Stirrer 220V includes 2 cord sets:

- United Kingdom/Singapore cord set and
- European cord set

Description

- | | |
|-----------------------------|--------------------------|
| ① Ergonomic Carry Handle | ① Splash Guard |
| ② Rubber Gasket | ② Housing - Molded Front |
| ③ Button - Power ON / OFF | ③ Housing - Cast Back |
| ④ Button - Increase RPM | |
| ⑤ Button - Decrease RPM | |
| ⑥ Four Digit LED Display | |
| ⑦ 1/2" Diameter Support Rod | |
| ⑧ Chuck | |

Technical Data

Speed Range	40 - 2000 rpm
LED Display Accuracy	± 1 rpm or 3% of reading, 4 digit
Maximum Torque	70 N-cm, (100 in-oz)
Input to Power Supply	AC 100-240 volts, 50/60 Hz, 125-175 VA
Electrical Output from Power Supply to Stirrer	DC 24 V, 2.7 A
Mechanical Output Power	1/15 hp, 50 watts, (Brushless DC Motor)
Weight	Stirrer Power Supply
	3.0 kg (7 lbs) 0.5 kg (1 lb)
Maximum Volume	25 litres (6.6 US gallons, 5.3 Imp Gallons)
Maximum Viscosity	15,000 mPa.s (15000 cps) - e.g. sour cream
Molded Front Housing	Flame retardant, chemical resistant, high impact nylon. Recyclable (PA).
Cast Back Housing	Epoxy coated zinc
Combined Housing	Splash Proof
Key Pad	Waterproof, chemical resistant polyester
Power Cord Lengths	1 M (3.2') from Stirrer to Power Supply-DC 1 M (3.2') from Power Supply to wall-AC
Adjustable Chuck	Fits onto 13mm (33/64 inch) Ø output shaft of stirrer. Adjusts up to 9.5 mm (3/8 inch).

Technical specifications subject to change without notice

Installation In The Laboratory

Installing the support bar

Remove the allen screw and lock washer from the support bar. Insert the support bar into the back of the stirrer so that the holes line up for the allen screw. Secure with allen screw and lock washer and tighten with allen key.

Clamp the stirrer to a Fisher Scientific stand on a stable, level surface. See Safety Instructions.

Connect stirrer to supplied Fisher Scientific power supply. Connect power supply to a grounded electrical outlet.

Impeller attachment

Insert the impeller shaft into the open chuck and adjust the impeller to the desired height. The hollow shaft allows the impeller shaft to come through the top of the stirrer if necessary.

CAUTION: Operating a freely rotating impeller shaft through the top of the stirrer is not safe. Do not allow the impeller shaft to touch or extend through the rubber gasket while the motor is powered on. The gasket is for sealing purposes only. Tighten the chuck's outer ring by hand and then secure with the chuck key.

Environmental Operating Conditions

The stirrer must operate in the following conditions:

Indoors

1. Altitudes up to 2000 m (6500 ft.).
2. Temperatures from 5°C to 40°C (41°F to 104°F).
3. Maximum relative humidity 80% for temperatures up to 31°C (88°F) decreasing linearly to 50% relative humidity at 40°C (104°F).
4. Pollution degree 2 in accordance with IEC 664
“Normally only non conductive pollution occurs. Occasionally, however, temporary conductivity caused by condensation must be expected.”

General Operating Instructions

Make sure the stirrer is securely mounted as specified in the installation instructions. Ensure the mixing paddle or impeller is securely attached to the desired position.

The stirrer is now ready for use to mix liquids and liquid/solid solutions. Appropriate impeller and container specifications must be determined by the operator, ensuring that all safety instructions are followed.

Button Operating Functions

Power Button - Turns the stirrer motor on and off.

Increase Button - Increases the speed by one rpm.
If pressed within five seconds of the stirrer being turned on it will stir at the last speed before powered off.

Decrease Button - Decreases the speed by one rpm.
If pressed within five seconds of the stirrer being turned on it will stir at 40 rpm.

Safety is built into this stirrer. When it is turned on it pauses for 5 seconds before it starts to mix. The last speed used is displayed on the LED. The user has three options.

1. Do nothing. The stirrer will start to stir at 40 rpm after five seconds.
2. Press the  and the stirrer will start to stir at the last speed used before powered off. This speed is displayed.
3. Press the  and the stirrer will start to stir at 40 rpm.

LED Display Displays desired rotating speed in rpm. A “soft start” has been incorporated into the stirrer. When the increase or decrease buttons are pressed and held, the speed ramps smoothly so that getting to the speed you require is quick and easy.

Display Codes

OL Indicates an overload condition. The load is too great for the stirrer and must be reduced. To reset stirrer turn power off and back on.

OH Indicates an overheat condition. The operating conditions have caused the stirrer to overheat. Reduce stirrer loads and / or environmental temperature to continue mixing. To reset stirrer turn power off, allow stirrer to cool and then turn back on.

Routine Cleaning, Service, Transport & Storage

Cleaning

Using a mild soap and water solution, wipe the housing and control panel of the stirrer with a damp cloth.

Note :

Do not use chlorine bleach, chlorine-based cleanser, abrasives, ammonia, steel wool or scouring pads with metal content or similar harsh solvents or abrasives. These may damage the surface of the instrument.

Service

Do not attempt to service or repair a Fisher Scientific stirrer. **If the stirrer housing is opened the warranty becomes void.** Contact Fisher Scientific for return authorization and return instructions.

Transport and Storage

1. Clean stirrer so that it is free from any materials which may be harmful to the health. Provide a MSDS or other safety data sheets for all chemicals used in conjunction with this stirrer where appropriate.
2. Place the unit and its parts into the original packing or a container with the necessary protection to prevent damage during transport. Seal the original packing or container with packing tape.
3. Store the packed unit in a dry place.

CAUTION :

Failure to clean, service, and handle the stirrer as outlined can cause it to fail which could impair the safety of the user.

Warranty and Product Liability

Warranty

Fisher Scientific warrants the herein described product of its manufacture for three years from receipt of the equipment by the purchaser, against defects in material and workmanship.

This limited warranty covers parts and labor.


If the stirrer housing is opened the warranty becomes void.

This warranty may only be altered by a specifically published amendment. No individual has authorization to alter the provisions of this warranty policy or its amendments. Expendable items are not covered by this warranty. Damages due to corrosion or accidental breakage are not covered.

Product Liability

Under no circumstances shall Fisher Scientific be liable for indirect, consequential or special damages of any kind.

Fisher Scientific reserves the right to make technical changes without prior notice.

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© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.