# Instruction Manual

## Fisher Scientific Power Overhead Stirrer

English  
page 3 to 16

Français  
page 17 to 30


Español  
page 31 to 44

Deutsch  
page 45 to 58

FSR3000MN002  
Rev00

ENGLISH . . . . .	3
<i>Safety instructions</i> . . . . .	4
<i>General instructions</i> . . . . .	5
<i>Stirrer approvals</i> . . . . .	6
<i>Packing list</i> . . . . .	6
<i>Description</i> . . . . .	7
<i>Technical data.</i> . . . . .	8
<i>Torque capabilities</i> . . . . .	9
<i>Specifications</i> . . . . .	10
<i>Keypad membrane</i> . . . . .	10
<i>Installation in the laboratory.</i> . . . . .	11
<i>Environmental operating conditions.</i> . . . . .	11
<i>General operating instructions.</i> . . . . .	12
<i>Zero torque feature</i> . . . . .	14
<i>Cleaninh and service</i> . . . . .	15
<i>Warranty and product liability.</i> . . . . .	16

## Safety Instructions

1. All operators must be familiar with the stirrer and should read this entire manual.

2. The stirrer must be securely fixed to a stable support. Mount to a Fisher Scientific Stand (Part Number 14-500-354 / 1298-7096) with a Fisher Scientific heavy duty clamp (Part Number 14-500-357 / 1290-7106). If other stands or clamps are used, the unit must be secure and can not fall if tipped up to ten degrees from the vertical.

3. The stirrer must be supplied with the rated voltage. See serial plate.

4. CAUTION: This is not an explosion proof stirrer. Do not use with highly flammable or explosive materials.

5. Spinning paddles or impellers can cause severe personal injuries. Operators must use extreme care and good judgement when mixing at any speed. The stirrer has higher torque capabilities than conventional stirrers.

6. All mixing paddles and impellers must be in good condition with straight shafts. If the stirrer vibrates at high speeds check the paddle shaft for damage and repair or replace it.

7. Extreme care must be taken when mixing chemicals so that no chemicals are splashed outside the mixing vessel. Care must be taken when changing to faster mixing speeds. On power up, the unit will display its set speed and will climb to that speed when the start/pause button is pushed. Always start at lowest speed if unsure of maximum safe speed.

8. Ensure that the mixing impeller does not contact the containing vessel.

9. Do not operate while standing in water. Keep the unit dry and do not immerse any part, except the mixing paddle, into any liquids. Protect from splashing.

10. Ensure that no loose clothing, jewelry, or hair can become entangled in any rotating parts. A fast spinning chuck can cause injury to an operator. Use chuck guard provided when stirrer is in use.

11. Power can be interrupted to the stirrer by pressing the power button or by disconnecting the mains cord. If rotating, this will cause the stirrer to stop and will disconnect power to the stirrer's internal drive circuit.

12. Shaft rotation can also be stopped by pressing the start/pause button. This does not disconnect power to the stirrer's internal drive circuit.

13. Wear safety goggles and suitable clothing when operating the stirrer.

14. Repairs must be carried out by authorized technicians only.

## General Instructions

1. Read all safety instructions.
2. Unpack the stirrer carefully. Ensure it is not damaged and verify the contents against the packing list. If the stirrer is damaged or the contents do not correspond with the packing list, notify your supplier immediately.
3. Read the manual carefully! Make sure that all users read this manual thoroughly before operating the stirrer.
4. This instruction manual should always be kept readily available for reference.
5. For further questions or inquiries, please contact:

### **Fisher Scientific**

United States

Customer Service: 1-800-766-7000

To Fax an Order: 1-800-926-1166

To Order online: [www.fishersci.com](http://www.fishersci.com)

Canada

Customer Service: 1-800-234-7437

To Fax an Order: 1-800-463-2996

To Order online: [www.fishersci.ca](http://www.fishersci.ca)

Europe

Austria: +43(0)800-20 88 40

Belgium: +32 (0)56 260 260

Denmark: +45 70 27 99 20

Germany: +49 (0)180 5258221

Ireland: +353 (0)1 885 5854

Italy: +39 02 950 59 478

Finland: +358 (0)9 8027 6280

France: +33 (0)3 88 67 14 14

Netherlands: +31 (0)20 487 70 00

Norway: +47 22 95 59 59

Portugal: +351 21 425 33 50

Spain: +34 902 239 303

Sweden: +46 31-68 94 30

Switzerland: +41 (0)56 618 41 11

UK: +44 (0)1509 555 500

# Stirrer Approvals

North America	 <b>CAN / CSA 61010-1 &amp; 61010-2-51</b> <b>IEC 61010-1 &amp; 61010-2-51</b>
Europe	 <b>CE, IEC 61010-1 &amp; 61010-2-51,</b> <b>EMC EN50081-1, EN50082-1</b> <b>(EN61000-4-2,</b> <b>EN61000-4-3, EN61000-4-4,</b> <b>EN61000-4-5)</b> <b>RoHS</b> <b>WEEE</b>

CSA Canadian Standards Association  
 IEC Low Voltage Act  
 EMC Electromagnetic Compatibility Generic Immunity Standard  
 RoHS Restriction of Hazardous Substances (EU)  
 WEEE Waste Electrical and Electronic Equipment (EU)


## Packing List

Compare the contents of the shipment to the items on the following list to ensure all parts are received with the unit. Do not discard the container and packing materials until all parts are accounted for. It is recommended to keep the original packaging for transportation and storage.

- Power Overhead Stirrer (115V Cat No. 14-500-210 or 220V Cat No. 1529-3797 )
- Chuck with key
- Chuck Guard
- Calibration Certificate
- Instruction Manual

Includes:

- United Kingdom/Singapore cord set and
- European cord set

## Description

This stirrer is specifically designed for professional use in the laboratory and for use in industry to stir numerous liquids with various properties and applications.

The stirrer rotates at 20 - 3000 rpm with a maximum torque of 341 Ncm (30 in-lbs)

The stirrer uses a digitally controlled high efficiency brushless DC motor which yields high torque over a very wide speed range. The speed is kept constant with continuous feedback to the motor. The stirrer is designed for reliability and durability. There are no wear parts that require user service during the normal life of the stirrer. There are no user serviceable parts.

The stirrer incorporates smart logic for automatic transmission control. This logic determines the appropriate speed range at which the stirrer should operate. See Fig 1 and Fig 2 under Torque Capabilities. The automatic selection can be overridden by pressing the speed button.

With higher solution viscosities, more torque is required to stir at the same speed. When the torque demands are outside the instrument's range, the control logic slows down the stirrer until the torque comes within the instrument's range. If the torque demand is too great at the stirrer's lowest possible speed, the unit will automatically shut down.

The stirrer is designed for safety. The unit comes on in "pause" mode once power is activated. After a power interruption the stirrer remains off. The stirrer's speed will gradually ramp up to the desired set point after the start/pause button is pressed.

The stirrer has integrated electronic overload protection which ensures no damage can occur to the electronics, motor or transmission when an overload situation occurs, such as too great a torque demand or unit overheating.

A chuck guard is supplied to place over the chuck. It protects the user from the spinning chuck and also helps to protect the chuck from potential splashes.

# Technical Data

Voltage	120 volts @ 60 Hz	220 volts @ 50/60 Hz
Installation Category (in accordance with IEC664)	II	II
Maximum Current Draw	3 Amps	3 Amps
Maximum Output Power	1/5 hp, 150 watts	
Speed Range	20-3000 rpm	
Speed Readout	LCD 4 digit display	
Speed Readout Accuracy	+/- 1 rpm or +/- 1% of reading	
Torque Readout Accuracy	+/- 12 N-cm (1 in-lb) or +/- 5% of reading	
Maximum Viscosity	50,000 mPas (50,000 cPs)	
Maximum Volume	60 litre (15.85 US Gal.)	
Maximum Operating Noise	70 dB(A) at 1 meter	
Motor Type	Brushless DC	
Drive	2 stage transmission	
Hollow drive shaft	Nitrated P .30 with 10 mm (0.39") diameter bore for stirring impellers 13 mm (0.51") outside diameter for drive chuck	
Chuck Guard	Length 67 mm, (2.64") Diameter 54 mm, (2.13") Molded of clear chemical resistant polyurethane.	
Seals	Butyl Rubber IP 42	
Holding chuck	3 prong stainless steel. Maximum opening 10 mm (0.39") diameter	
Mounting Rod Length	140 mm, (5 1/4")	
Stirrer Weight	5 kg, (11 lbs.)	

Protection measures:		0
Class	Protection device	Recommended use
0	No additional protection	Used only on equipment where failure of the controls does not represent a hazard.
The operating conditions should be monitored at regular intervals.		

# Torque Capabilities

Fisher Scientific	Low Speed				High Speed			
	Min Speed rpm	Max Speed rpm	Peak Torque Ncm	Cont. Torque* Ncm	Min Speed rpm	Max Speed rpm	Peak Torque Ncm	Cont. Torque* Ncm
Power Overhead Stirrer	20	584	341	190	206	3000	68	43

\*Continuous torque at 50% of rpm range and 25 degrees Celcius.

Torque Conversion Table		
In-lb to In-oz	Multiply by 16	50 in-lb = 800 in-oz = 568 N-cm
In-lb to N-cm	Multiply by 11.36	30 in-lb = 400 in-oz = 341 N-cm
N-cm to In-oz	Multiply by 1.416	15 in-lb = 240 in-oz = 170 N-cm

Fig 1




Fig 2


# Specifications

Fig 3

# Keypad Membrane

Fig 4

## Installation in the Laboratory

1. The stirrer must be securely fixed to a stable support - e.g. clamped to an original Fisher Scientific stand that is on a stable, level surface. See Safety Instructions.

CAUTION:

This is not an explosion-proof stirrer. Do not use with highly flammable or explosive materials.

WARNING:

Be sure that the unit's voltage is the same as the supplied voltage. Information on the unit's voltage and current are located on the serial plate.

This unit must always be connected to a grounded electrical outlet.

2. Chuck removal and attachment:

Loosen the set screw with the allen key found on the chuck key. Remove chuck. To attach align the set screw with the pilot hole on the shaft. Tighten the set screw with the allen key found on the chuck key. Make sure that the set screw is firmly tightened. Open the chuck fully to prepare for impeller attachment.

3. Impeller attachment:

Insert the impeller shaft into the open chuck and adjust the impeller to a maximum length of 25 cm (10 inches) [of a 40 cm (16 inch) shaft] to reduce vibration. The hollow shaft allows the impeller shaft to come through the top of the stirrer if necessary. Tighten the chuck's outer ring by hand and then secure with the chuck key. Install the chuck guard over the chuck.

## Environmental Operating Conditions

The stirrer must operate in the following conditions:

### Indoors

1. Altitudes below 2000 meters.
2. Temperatures from 5 °C to 40 °C
3. Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C.
4. Mains supply as indicated under the technical specifications.
5. Pollution degree 2 in accordance with IEC 664  
"Normally only non conductive pollution occurs.  
Occasionally, however, temporary conductivity caused by condensation must be expected"

# General Operating Instructions

-  1. Make sure that the stirrer is securely mounted as specified in the installation instructions. Also make sure that the mixing paddle or impeller is securely attached in the desired position.
-  2. Connect the power cord to a grounded supply.
-  3. The stirrer is now ready for use to mix liquids and liquid/solid solutions. The operator, ensuring that all safety instructions are followed must determine appropriate impeller and container specifications.
-  4. The stirrer can operate continuously if the torque demand falls below the continuous torque operation line (see Torque Capabilities Fig 1 and Fig 2). Should the ambient temperature be above 25 °C, then the 40 °C continuous operation line must be used.
-  5. If the stirrer is used for an extended period of time above the appropriate continuous temperature operation line (see Torque Capabilities Fig 1 and Fig 2), then it may shut down. The stirrer may be started after it has cooled down.

## Button Operating Functions

- **Power** - when power button is pressed for the first time stirrer is activated in “pause” mode. The set speed is displayed. RPM LED will be flashing indicating this is a set speed.
- When stirrer is on the Power LED is on.
- When the **Power** button is pressed a second time, the stirrer is turned off. The last set speed is remembered.
- **Start/Pause** - stops the stirrer when operating.
- **Start/Pause** returns the stirrer to the previous condition and previous display or starts the stirrer if stopped.
- Pause LED is on when stirrer is paused and set speed is displayed.
- **Display** - When running, toggles display between speed (RPM) and torque (N-cm). The LEDs indicate which parameter is displayed.
- When paused, the display will always show set speed.
- **Speed** - determines speed range (high speed/low speed) and overrides automatic choice when allowable. The LEDs indicate which speed range is operating.
- +/- (**RPM**) Buttons - Press buttons once to increase / decrease set speed by one rpm. Press and hold to increase / decrease speed at a greater rate.
- When unit is in pause, user can change the set speed by pressing the buttons.
- RPM LED will flash when displaying desired set speed and will be constant when displaying actual speed.
- As unit is auto ramping to a faster speed the user can stop the increase in speed by pressing the down arrow button. This changes the set speed to the present actual value.
- As unit is auto ramping to a slower speed the user can stop the decrease in speed by pressing the up arrow button. This changes the set speed to the present actual value.

## Indicator Lights

- High Speed / Low Speed LED's - These indicate the stirrer's set range.
- Overload LED - This indicator will flash when an overload condition is approached. This indicator will remain on when an overload condition is reached and the stirrer is shut down.

Overload conditions will occur when:

- The torque demand is outside the unit's capabilities (solution is too viscous).
- The unit overheats due to extended operation at very high torque and/or speed demands.

When the overload condition is due to torque demand, the stirrer will automatically reduce speed to reduce the load. If reduction in speed is not sufficient to keep the load in the stirrer's capable range, then the stirrer will shut down.

When the overload condition is due to overheating, the stirrer will not automatically reduce speed. If the overload light is flashing and the speed has not reduced from the set point, the load should be reduced by the operator.

In either case, if the stirrer completely shuts down due to an overload condition, it can be reset by powering off the stirrer. If the stirrer has overheated then it must cool down before it can be reset.

# Zero Torque Feature

## Theory

Torque (turning force) is a measure of the force on an object such as a shaft or a propeller. For example, pushing or pulling the handle of a wrench connected to a nut or bolt produces a torque (turning force) that loosens or tightens the nut or bolt. The units used to quantify torque are Newton-meters, Newton-centimeters, or foot pounds.

In **Mixing Applications** the magnitude of torque depends on three quantities: First, the speed in revolutions per minute; second, the diameter of the propeller; and third, the viscosity of the material being mixed. Additionally, a small amount of torque is required to move the mixing impeller. In small bench top applications this additional or overhead torque is negligible and is generally independent of the force required to mix material.

## Zero Torque Feature

The zero torque feature allows the operator to null out the additional overhead torque of the stirring apparatus similar to a tare feature on a scale. The feature can also be used during the mixing process and allows the operator to see change in torque caused by the change in material viscosity.

## Zero Torque Operating Instructions

Run the stirrer at low speed for a ½ hour before zeroing torque. As the motor warms up the torque can change slightly. Setting the torque when the stirrer is warm will increase the accuracy.

1. Follow installation instructions. Before attaching the mixing impeller, press the power button, ramp to the desired speed and press the start/pause button to start the stirrer. Speed in rpm is displayed.
2. Press the display button. The stirrer displays actual torque. Note this cannot be done unless the stirrer is running.
3. Press the Zero Torque button. The display now reads zero and the N-cm LED will flash to indicate that the torque has been set to zero. Pause mixing. The RPM LED will flash to indicate the stirrer has been paused.
4. Place the impeller into the chuck and then into the mixing vessel and press the start/pause button to start mixing. The torque of the material being mixed will now be displayed.

Follow all safety instructions when mixing.

## Note:

At any time while mixing the torque display can be set to zero. Viscosity changes can then be monitored as solutions change during the mixing process. A negative number indicates a reduction in viscosity and a positive number indicates an increase in viscosity.

### Other button functionality

- If the stirrer is powered off when in zero torque mode (N-cm LED flashing), the stirrer will be in this mode when powered back on. If the stirrer is unplugged or power is interrupted, the stirrer is reset and the feature will not be active.
- The speed can be changed any time – while the stirrer is running or in pause. When the speed is changed, the speed in rpm is displayed while the speed is adjusted. If torque was displayed prior to the change in speed, the display will change back to displaying torque once the new speed is set.
- The zero torque feature can only be activated when the stirrer is running.
- When the zero torque button is pressed while the stirrer is running and rpm is displayed, actual torque will be displayed. Press the zero torque button again, the N-cm LED flashes indicating the stirrer is in zero torque mode.

## **Cleaning and Service**

### Cleaning

Using a mild soap and water solution, wipe the housing and control panel of the stirrer with a damp cloth. Do not use chlorine bleach, chlorine-based cleanser, abrasives, ammonia, steel wool, or scouring pads with metal content or similar harsh solvents or abrasives. These may damage the surface of the stirrer.

### Service

Do not attempt to service or repair a Fisher Scientific stirrer. **If the stirrer housing is opened, the warranty becomes void.** Contact Fisher Scientific for return authorization and return instructions.

Clean the stirrer so that it is free from any materials which may be harmful. Provide an MSDS or other safety data sheets for all chemicals used in conjunction with this stirrer where appropriate. Place the stirrer and its parts into the original packaging to protect and prevent damage during transport. Follow shipping directions provided by Fisher Scientific when return is authorized and RMA issued.

 **CAUTION:** Failure to clean, service, and handle the stirrer as outlined above can cause it to fail which could impair the safety of the user.

# Warranty and Product Liability

## Warranty

Fisher Scientific warrants the herein described product of its manufacture for three years from receipt of the equipment by the purchaser, against defects in material and workmanship.


This limited warranty covers parts and labor. In the event of a warranty claim, contact Fisher Scientific. If the cause is determined to be a manufacturing fault, Fisher Scientific will repair or replace all defective parts to restore the unit to operation.

This warranty may only be altered by a specifically published amendment. No individual has authorization to alter the provisions of this warranty policy or its amendments. Expendable items are not covered by this warranty. Damages due to corrosion, accidental breakage or submersion are not covered.

## Product Liability

Under no circumstances shall Fisher Scientific be liable for indirect, consequential or special damages of any kind.

Fisher Scientific reserves the right to make technical changes without prior notice.

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© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).