

# Laboratory identification Labels



Products that perform. Quality that lasts



Extensive material knowledge combined with chemical expertise enables Brady to develop durable labels that perform throughout sample storage and processing.

- B-7425 An all-round laboratory label: Track almost any sample in the lab. It is designed to stay attached and remain legible in a wide range of sample collection, sample processing, microbiology and cryogenic applications
- Brady's FreezerBondz™ B-492 Labels are designed specifically to identify laboratory samples for cryo storage. FreezerBondz labels can be applied to room temperature or to already frozen samples, albeit with a small overlap. The label will remain legible and will stick to the sample container during processing, even when exposed to extreme heat or chemicals like DMSO, Xylene and Ethanol
- StainerBondz™ Labels B-481 are designed specifically to identify slides and to resist exposure to harsh chemicals during the slidestaining process. They perform well in the hematoxylin and eosin slidestaining process and offer superior chemical resistance also useful in other laboratory contexts. StainerBondz™ can resist temperatures ranging from -80°C (freezer) up to 121°C (autoclave).

Cat. No.	Description	Pack Qty
15318982	Brady B-7425, BMP21 printer : M21-750-7425	1 Cartridge
17150191	Brady B-7425, BMP51 printer : M-31-7425	250 Labels/Pk
17140191	Brady B-7425, M611 printer : M61-98-7425	500 Labels/Roll
17150211	Brady B-7425, J2000 printer : J20-88-7425J	2000 Labels/Roll
17130181	Brady B-7425, i3300 printer : B33-88-7425	2000 Labels/Roll
17180161	Brady B-7425, BBP12 printer : THT-152-7425-2-SC	2000 Labels/Roll
12662996	FreezerBondz™ B-492, BMP51 printer : M-131-492	180 Labels/Pk
17120421	FreezerBondz™ B-492, M611 printer : M61-19-492	250 Labels/Roll
15354158	FreezerBondz™ B-492, i3300 printer : B33-59-492	2500 Labels/Roll
15529275	FreezerBondz™ B-492, BBP12 printer : THT-250-492-1.5-SC	1500 Labels/Roll
15490945	StainerBondz™ Labels B-481, BMP51 printer : M-194-481	320 Labels/Pk
17150411	StainerBondz™ Labels B-481, M611 printer : M61-195-481	500 Labels/Roll
13573330	StainerBondz™ Labels B-481, i3300 printer : B33-194-481	2500 Labels/Roll
12702447	StainerBondz™ Labels B-481, BBP12 printer : BPT-194-481	2500 Labels/Roll

[Visit our brand page](#)

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk



© 2022 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.