

# Lysing Beads for Homogenisers

Fisherbrand offers a selection of lysing beads and tubes for processing a wide range of samples.

A variety of bead media are available for different sample types and bead tubes are offered pre-filled. Fisherbrand bulk beads and empty bulk tubes are also available. Although best results are achieved with the Fisherbrand Bead Mills 4 & 24, Fisherbrand bead tubes are also compatible with most other commercially available bead mills.

## Fisherbrand Grinding Beads Selection Guide

Sample Type	Bead Media			
	Glass	Ceramic		Metal
	0.5mm	1.4mm	2.8mm	2.4mm
Microbes and spores	•			
Soft animal tissues (brain, liver, kidney)		•	•	
Fibrous animal tissues (heart, muscle, skin)			•	
Plant materials			•	
Difficult biological samples (hair, bone, seeds etc)				•

## Fisherbrand Prefilled Bead Tubes

Cat. No	Bead media	Bead Size	Tube Size	Pack qty
15545799	Glass	0.5mm	2mL	50
15555799		1.4mm	2mL	50
15565799	Ceramic	2.8mm	2mL	50
15585799		1.4mm	7mL	50
15595799		2.8mm	7mL	50
15535799	Metal	2.4mm	2mL	50
15575799		2.4mm	7mL	50



## Fisherbrand Bulk Beads

Cat. No	Bead media	Bead Size	Tube Size	Pack qty
15515809	Ceramic	1.4mm	2mL	325g
15525809		2.8mm	2mL	325g
15505809	Metal	2.4mm	2mL	500g



## Fisherbrand Empty Bulk Tubes

Cat. No	Description	Tube Size	Pack qty
15535809	Reinforced bulk tubes (supplied empty without beads) with screw caps and O-rings	2mL	500
15545809	Non-reinforced bulk tubes (supplied empty without beads) with screw caps and O-rings	2mL	500



Find out more at [eu.fishersci.com](http://eu.fishersci.com)

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Czech Republic: +420 466 798230  
 Denmark: +45 70 27 99 20 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14  
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478  
 Lithuania: +370 5 244 4442 Netherlands: +31 (0)20 487 70 00 Norway: +47 22 95 59 59  
 Portugal: +351 21 425 33 50 Spain: +34 902 239 303 Sweden: +46 31 352 32 00  
 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

Other Central & Eastern Europe: +49 2304 932 715 Middle East: +44 1509 555939  
 Africa (Anglophone): +44 1509 555939 Afrique (Francophone): +33 3 8867 5327

© 2019 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at [fishersci.com/trademarks](http://fishersci.com/trademarks).